

# Sensors and Digitization

## Moving Object Imaging

### Lab 1 Report

### GrTP1A

Daniel Sileshi Asfaw  
Emre Ozan Alkan

4 December 2013

## 1 Introduction

### 1.1 Objective

The main goal was to understand how a simple camera can be transformed into a polarization state measurement system. This lab introduces two different systems:

1. A simplified polarization imaging set-up that consist of a manually rotating polarizer placed in front of the sensor.
2. A contrast polarization measurement system that uses a Twisted Nematic Liquid crystal.

### 1.2 Equipment

In the lab room, we had following equipment to use:

- PC Computer
- Frame Grabber - IEEE 1394
- Camera - Allied Vision Technologies GUPPY + one lens + Video Cable
- Arcoptix switchable polarization rotator 0-90 degree (Twisted Nematic Liquid Crystal) + Arcoptix USB LC Driver
- Two linear polarizers
- Four mounting posts and four post holders
- Lighting Device + Polarized Ring

### 1.3 Software

We used the following software in the lab computers:

- National Instruments "Measurement & Automation Explorer" Software
- Arcoptix USB LC Software
- Matlab or National Instruments LabVIEW

In this laboratory, we performed three different experiments; the first experiment was grabbing images using one linear polarizer, the second was contrast polarization measurement using two linear polarizer also we have experimented using liquid crystal polarizer on the middle of the two linear polarizer, the last experiment was to study one application of polarization, diffuse specular reflection.

## 2 Simplified Polarization Imaging

### 2.1 Wolffs Method




Figure 1: First Setup

We first setup all these imaging equipment as in figure [1]. We adjusted the camera position will capture computer screen as a part of the grabbing scene. In order to avoid reflection seen from appearing on the graded image we position the camera as close as possible to the polarizer. The "Measurement & Automation Explorer" Software is opened to start grabbing an image. We start grabbing an image by changing the angle of the polarizer. From the grabbed image we observe that the computer screen's contrast value was caning while we were changing the angle of the polarizer. The reason of this phenomenon is that computer screen is already polarized.

Here is the image grabbed with different polarizer orientations(0, 45 and 90 degrees)




Figure 2: 0 - 45 - 90

We used MATLAB to calculate the special case of least mean square method (Wolff's method) to find the following parameters through three grabbed images as in equation below:

$$\begin{cases} I &= I_0 + I_{90} \\ \tan 2\varphi &= \frac{I - 2I_{45}}{I - 2I_{90}} \\ \rho &= \frac{\sqrt{(I - 2I_{45})^2 + (I - 2I_{90})^2}}{I} \end{cases}$$

Figure 3: Wolff's Method

Those three values are HSV color image representation of the image. Where H is angle of polarization S is degree of polarization and V is the intensity.




Figure 4: Wolff's Method Gray Scale Results






Figure 5: Wolff's Method HSV Result

## 2.2 Least Mean Square Method

We took 8 photos with degrees 0, 45, 90, 135, 180, 225, 270, 315, respectively:



Then we applied Least Mean Square Method for images as we seen in class. Here is the results:




Figure 6: Least Mean Square Method Result

And here is the sinusoidal relationship of the images:




Figure 7: Least Mean Square Method Sinusoidal Relationship

### 3 Contrast polarization measurement

#### 3.1 2 Linear Polarizers




Figure 8: 2 Linear Polarizer Setup

We setup all the imaging equipment as in figure above. We start by adjusting polarization angle to of the two linear polarizer to zero. We observed that output image is darker but more readable than one

from first experiment at same orientation using a single polarizer. Also intensity decreased a lot because we using 2 polarizer as a result of malus' law.

### 3.2 Introducing Liquid Crystal




Figure 9: Introducing Liquid Crystal

The liquid crystal is put between the polarizer. This switchable polarization rotator works by controlling the amount input voltage using Arcoptix software. It works in binary mode. It rotates to 90 degree when input voltage is applied and zero degree when input voltage is 0V. This set up helps to automate the rotation when experimenting effect of polarization using two linear polarizers.

### 3.3 Diffuse Specular Reflection




Figure 10: Using Polarized Lighting Source

We setup all the imaging equipment as in figure above by putting Polarized Ring with Lighting Device in front of the polarizers. And we started snapping two different images of two different objects (hand and metal object) at two different polarization angles (0 and 90 degrees).




Figure 11: Removing Specular Reflection

We observed that specular reflection (smooth or shiny surface of an object) can be removed by setting the polarization rotator to 0V (shiny metal example), but diffuse reflection like hand, can't be removed at any setting.

## 4 MATLAB CODES

### 4.1 Wolffs Method

Listing 1: Wolffs

```

1 image_zero = imread('0.png');
2 I_0 = double(image_zero);
3
4
5 image_45 = imread('45.png');
6 I_45 = double(image_45);
7
8 image_90 = imread('90.png');
9 I_90 = double(image_90);
10 %===== Wolff's method =====
11 %intiallize the out put images to zero first
12 I = zeros(size(I_0));
13 Phi = zeros(size(I_0));
14 Ro = zeros(size(I_0));
15
16 for i=1:size(I_0,1)
17     for j=1:size(I_0,2)
18         for k=1:size(I_0,3)
19             I(i,j,k) = I_0(i,j,k) + I_90(i,j,k);
20             Phi(i,j,k) = (((I(i,j,k)-2*I_45(i,j,k))/(I(i,j,k)-2*I_90(i,j,k))));
21             Ro(i,j,k) = (sqrt(((I(i,j,k)-2*I_45(i,j,k)).^2) + ((I(i,j,k)-2*I_90(i,j,k)).^2))/I(i,j,k));
22         end
23     end
24 end
25
26
27
28 %create the 3D color image by concating the above three results and
29 %display the image
30 Image_out = uint8(cat(3,atan(Phi),Ro,I));%this image is hsv image
31
32 figure,subplot(131), imshow(image_zero);title('Image at zero degree of pol');
```

```

33 subplot(132),imshow(image_45);title('Image at 45 degree of pol');
34 subplot(133),imshow(image_90);title('Image at 90 degree of pol');
35
36 figure (2),subplot(131), imshow(uint8(I));title('Intensity');
37 subplot(132),imshow(uint8((Phi)));title('Angle of polarization');
38 subplot(133),imshow(uint8(Ro*100));title('Degree of Polarization');
39 figure(3), imshow(Image_out);colormap(hsv);title('out put image color image');

```

## 4.2 Least Mean Square Method

Listing 2: Least Mean Square Method

```

1 Im(1, :, :) = double(imread('0.png'));
2 Im(2, :, :) = double(imread('45.png'));
3 Im(3, :, :) = double(imread('90.png'));
4 Im(4, :, :) = double(imread('135.png'));
5 Im(5, :, :) = double(imread('180.png'));
6 Im(6, :, :) = double(imread('225.png'));
7 Im(7, :, :) = double(imread('270.png'));
8 Im(8, :, :) = double(imread('315.png'));
9
10 image_zero = imread('0.png');
11 Im_0 = double(image_zero);
12
13 [m n] = size(Im_0);
14
15
16 Im_out = zeros(m,n,3);
17 for i=1:m
18 for j=1:n
19 A = [];
20 Y = [];
21 for k=1:8
22
23 A = [ A; 1 cos(2*k*45) sin(2*k*45) ];
24 Y = [ Y; Im(k,i,j) ];
25
26 end
27
28 X = 2 .* (inv(A' * A)) * (A' * Y);
29 I = X(1);
30 Ro = sqrt(X(2).^2 + X(3).^2) ./ X(1);
31 Phi = atan(X(3) ./ X(2));
32 Im_out(i,j,:) = [ I Ro Phi];
33 end
34 end
35
36 figure(6);subplot(131); subimage(Im_out(:,:,2)); title('degree of polarization');
37 subplot(132); subimage(Im_out(:,:,3)); title('angle of polarization');
38 subplot(133); subimage(uint8(Im_out(:,:,1))); title('intensity');
39
40 % Sinusoidal Relationship
41 yaxis=[Im_0(162,520) Im_45(162,520) Im_90(162,520) Im_135(162,520) Im_180(162,520) Im_225
(162,520) Im_270(162,520) Im_315(162,520)];
42 xaxis = [0 45 90 135 180 225 270 315];
43 figure(5);plot(xaxis,yaxis,'—rs','LineWidth',2,...
44 'MarkerEdgeColor','k',...
45 'MarkerFaceColor','g',...
46 'MarkerSize',10);
47 title('sinusoidal relationship');
48 xlabel('Angles');
49 ylabel('ROI');

```



### 4.3 Diffuse Specular Reflection

Listing 3: Specular Reflection Removal

```
1 % Diffuse Specular reflection
2 parell = imread('Metal_0V.png');
3 Orto = imread('Metal_1V.png');
4
5 Total_Intensity = parell + Orto;
6 Contrast = parell - Orto;
7 Contrast_ratio = Contrast ./ Total_Intensity;
8
9
10 subplot(1, 3, 1), imshow(uint8(Total_Intensity)); colormap(gray); title('Total Intensity')
11 ;
12 subplot(1, 3, 2), imshow(uint8(Contrast)); colormap(gray); title('Image Contrast');
13 subplot(1, 3, 3), imshow(uint8(Contrast_ratio), []); colormap(gray); title('Contrast ration
    ');
```